

INSTRUCTIONS FOR USE

DEVICE SYSTEM NAME:

GREENS SURGICALS (P) LTD NON-STERILE INTRAMEDULLARY NAILING SYSTEM

Indications for Use: **Tibial Nail**

GREENS SURGICALS (P) LTD. Tibial Nail System is intended to stabilize fractures of the proximal and distal tibia and the tibial shaft; open and closed tibial shaft fractures; certain pre- and post-isthmic fractures; and tibial malunions and non-unions.

The tibial nails can be of Interlocking tibial nail, cannulated tibial nail and distal tip distal locking tibial nail types.

Indications for Use: **Unreamed Femoral Nails**

GREENS SURGICALS (P) LTD. Retrograde/Antegrade Femoral Nail System is intended to stabilize fractures of the distal femur and the femoral shaft, including supracondylar fractures, including those with intra-articular extension; ipsilateral hip/shaft fractures; ipsilateral femur/tibia fractures; femoral fractures in multiple trauma patients; fractures proximal to total knee arthroplasty; fractures distal to a hip implant; fractures in the morbidly obese; fractures in osteoporotic bone, impending pathologic fractures; and malunions and nonunions.

Indications for Use: **Supracondylar Nails**

The GREENS SURGICALS (P) LTD. supracondylar intramedullary nails are intended for use in fixation of various types of fractures of supracondylar or distal region of the femur as well as for proximal fractures in conjunction with supracondylar fractures.

Indications for Use: **Proximal Femoral Nail, Proximal Femoral Nail-Long**

The GREENS SURGICALS (P) LTD. proximal femoral nail is intended to treat stable and unstable proximal femoral fractures including pertrochanteric fractures, intertrochanteric fractures, high subtrochanteric fractures.

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 1 OF 10

INSTRUCTIONS FOR USE

Indications for Use: **Elastic Nails**

The GREENS SURGICALS (P) LTD. Elastic Intramedullary Nail System is indicated for fixation of diaphyseal fractures where the canal is narrow or flexibility of the implant is important. This includes upper extremity fractures in all patients and lower extremity fractures in pediatric or small-stature patients. This system is also intended to treat metaphyseal and epiphyseal fractures, such as radial neck fractures and is intended for fixation of small long bones, such as carpal and tarsal bones.

MATERIALS: THESE DEVICES CAN BE MADE IN THE FOLLOWING MATERIAL GRADES

STAINLESS STEEL: 316L, 316LVM ALLOY

TITANIUM ALLOY in compliance with ASTM F 136.

CONTRA INDICATIONS:

The implant should not be used in a patient who has currently, or who has a history of:

- ▶ Local or Systemic acute or chronic inflammation.
- ▶ Active infection or inflammation.
- ▶ Suspected or documented metal allergy or intolerance.
- ▶ Symptomatic Arthritis.

Caution: "Federal law restricts this device to sale by or on the order of a physician"

WARNING:

Serious Post-operative complications may occur from the use of implant in a patient who-

- ▶ **Lacks good general physical condition**
- ▶ **Has severe osteoporosis**
- ▶ **Demonstrates anatomical or physiological anomalies.**
- ▶ **Has immunological response, sensitization or hypersensitivity for foreign materials.**
- ▶ **Do Not Re-Use. Do not re-process. The device is not designed for re-processing. Reprocessing of the medical device intended for single use devices may lead to degraded performance or a loss of functionality. The cleaning and disinfection / sterilization for re-processing of these single use devices has**

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 3 OF 10

INSTRUCTIONS FOR USE

Indications for Use: **Cannulated Humeral Nails**

The GREENS SURGICALS (P) LTD. Cannulated Humeral Nail System is intended to aid in the alignment and stabilization of humeral fractures to include:

- Diaphyseal fractures of the humeral shaft
- Fractures of the proximal humerus
- Proximal humeral fractures with diaphyseal extension
- Impending pathologic fractures
- Malunions and nonunions

Indications for Use: **Unreamed Humeral Nails**

The GREENS SURGICALS (P) LTD. unreamed humeral nails are indicated for the treatment of fractures of proximal humerus.

Indications for Use: **Universal Femoral Nails**

The GREENS SURGICALS (P) LTD. femoral nail is indicated for use in a variety of femoral fractures (Fig 1), such as:

- A. Comminuted fractures
- B. Segmental fractures
- C. Fractures with bone loss
- D. Proximal and distal fractures
- E. Nonunions
- F. Subtrochanteric fractures
- G. Intertrochanteric fractures

Indications for Use: **Recon Nails**

The GREENS SURGICALS (P) LTD. Recon Nail System is a fracture fixation device comprised of reconstruction nails and the related accessories such as washers, locking screws, set screws, end caps, and lag screws. The devices are intended to provide strong and stable internal fracture fixation with minimal soft tissue irritation of the tibia and femur. This device is utilized as an aid to healing, not as a substitute for normal intact bone and tissue

Indications for Use: **Ga-mma Nail**

The GREENS SURGICALS (P) LTD. Ga-mma Nail is indicated for use in fixation of femoral fractures occurring from the base of the femoral neck extending distally to a point proximal to intercondylar notch. Fracture types include basilar neck, intertrochanteric, subtrochanteric fractures and femoral shaft fractures.

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 2 OF 10

INSTRUCTIONS FOR USE

**not been validated nor any authentic information is available.
So re-process of the single use device is not allowed.**

Adequate knowledge of Surgical Techniques, proper selection and placement of implants and post-operative patient care are important considerations while performing an orthopaedic surgical procedure. Therefore it is essential that the implants are used by a qualified Medical Practitioner only.

NOTE: It is the responsibility of the surgeon to discuss, with the patient, the precautions, possible risks, warnings, consequences, complications and adverse reactions which may occur as a result of the surgical procedure and implantation of the device(s).

The patient should be informed that the life expectancy of the device is unpredictable and once implanted, successful results cannot be guaranteed if adequate care and precautions are not taken.

If adverse effects happen, it may necessitate re-operation, revision or removal surgery, arthrodesis or the involved joint and/or amputation of the limb

MRI COMPATIBILITY:

The GREENS SURGICALS (P) LTD. intramedullary nailing systems have not been evaluated for safety and compatibility in the MR environment. The GREENS SURGICALS (P) LTD. intramedullary nailing systems have not been tested for heating or migration in the MR environment.

I SELECTION OF THE PATIENT

During selection of patients the following factors have to be considered:

a) Immunological intolerance

Immunological intolerance may occur in some patients. Where material sensitivity is suspected, appropriate foreign body tests should be performed

b) Degenerative Diseases

In the case of patients suffering from degenerative diseases, this can get aggravated during implantation, and will decrease the expected life of the implant. In such cases surgery can be considered only as a temporary relief.

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 4 OF 10

INSTRUCTIONS FOR USE

c) Mental illness

Mental illness or schizophrenia may cause patients to ignore the limitations and precautions of the implanted material, leading to implants fracture and complication

d) Alcohol and Drug Addiction

The patients, addicted to alcohol and drugs, ignore, during the state of stupor or during the stage of withdrawal, the necessary precautions for the use of implants. This may result in complications of implants fracture.

e) Obesity

An obese patient produces abnormal stresses leading to increased load on the implant which ultimately results in failure of the implant.

f) Activity

If the patient indulges in an activity involving significant muscular strain in the implanted region, the result may be failure of the implant.

II SELECTION OF THE IMPLANT/IMPLANT SYSTEM

For successful implantation, the selection of the proper size, the shape and design of the implant are important.

The size and the shape of the human bones puts limitations of and strength of implants during fracture management and re-constructive surgery. The Doctor/Orthopedic Surgeon has to select the appropriate implant depending upon the application, which should also be free from any apparent corrosion or any manufacturing defect.

III HANDLING OF IMPLANT/IMPLANT SYSTEM

The GREENS SURGICALS (P) LTD. intramedullary nailing systems should conform to the shape of the bone. Any inequality should be adjusted by accurate shaping of the nail. To fit the bone is not recommended as the screws are instantly subjected to a strong expelling force.

Care should be taken that there are no scratches or distortions notches, sharp dents or reverse bends at the site of the screw hole.

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 5 OF 10

INSTRUCTIONS FOR USE

VIII CAUTION

Implant components from one manufacture should not be used with those of another. Implants from each manufacture may have metal and design differences so that the use in conjunction with different devices may lead to inadequate fixation or corrosion of the implant due to generation of peizo currents.

IX PACKAGING DISPOSAL

The packaging material of this device if swallowed, may cause choking Hazards. Therefore, it should be disposed of in such a way that it is out of reach of children and stray animals.

X STERILITY

GREENS SURGICALS (P) LTD. Intramedullary Nailing System are supplied Non-Sterile.

Check the integrity of the packaging and labelling before opening the pack.

Remove the device from the pack before sterilization.

Implants are recommended to be sterilized, using steam autoclaving process regularly used in the hospitals and clinics.

The following two methods have been validated and recommended by the company:

Method:	Steam Sterilization (Autoclaving)
Temperature	121 Degree Centigrade
Exposure Time:	15 Minutes

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 7 OF 10

INSTRUCTIONS FOR USE

These may cause defects on surface finish and result in improper bonding.

IV IMPORTANT

If the screw is not properly centred as it is being tightened, there will be a tendency for the countersink in the nail to force it to one side, thereby damaging the good threads already cut in the bone and producing a strain which will probably lead to necrosis and early loosening of the screw.

V IMPLANT REMOVAL

The surgeon must take the final decision on implant removal. It is recommended that an implant, used as an aid for healing should be removed once its services are over, particularly in younger and more active patients

VI POST OPERATIVE CARE

A patient must be made aware of the limitations of metallic implants and take precautions to avoid unnecessary stress to implant.

No partial weight bearing or non-weight bearing device can be expected to withstand the unsupported stresses of full weight bearing or excessive muscular activity when there is an un-united fracture. Bone union is a must and the patient must restrict his activities to assist in healing.

VII SECOND HAND IMPLANTS

Used implants, even if they appear un-damaged, may have internal and external Defects. It is possible that individual stress analysis each part fails to reveal the accumulated stress on the metals as a result of use within the body. This may ultimately lead to implant failure.

Every implant must be discarded after use and should never be re-used. It should be bent & then disposed of properly so that it becomes unfit for reuse. While disposing it of, it should be ensured that the discarded implant does not pose any threat to children, stray animals and environment.

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 6 OF 10

INSTRUCTIONS FOR USE

XI USE

Qualified & specialized trained surgeons only should use this device.






STORAGE

Store in a dry place.

PRODUCT INFORMATION DISCLOSURE:










GREENS SURGICALS (P) LTD. HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND MANUFACTURE OF THESE PRODUCTS. GREENS SURGICALS (P) LTD. EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GREENS SURGICALS (P) LTD. WILL NOT BE LIABLE FOR ANY INCIDENT OR CONSEQUENTIAL LOSS, DAMAGE OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THE DEVICES. GREENS SURGICALS (P) LTD. NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE DEVICES.

Details of various Symbols used in Labeling

	Do not reuse
	Batch code
	Date of manufacture YYYY-MM
	Consult Instructions for Use
	Manufacturer

SECTION : GREENS SURGICALS (P) LTD. /IFU/NAILS
REV.01 Dated: 17-07-2015 PAGE 8 OF 10

INSTRUCTIONS FOR USE

	Use by YYYY-MM
	Non-Sterile
	Catalogue number
	Keep Dry
	Keep Away from Sunlight
	Do not use if package is damaged
	EC Representative Address
	Material Grade Stainless Steel like 316L, 316LVM, Titanium TiAl6V4
	Quantity in the package

INSTRUCTIONS FOR USE

INFORMATION

Should any information regarding the implants or their use be required, please contact your representative or distributor or contact the manufacturers direct.
 Manufactured By:

Manufactured By:

GREENS SURGICALS (P) LTD.
 FDA: ISO 13485/2003 Certified
 Plot No: 500-610, GPOC, Seval Industrial Estate,
 Mariposa, Toboeka 391775 (Gujarat)
 Phone: +91-02667-264868, 264989, Toll Free: +91-02667-264860
 Web Site: www.greensurgicals.com
 E-mail: info@greensurgicals.com


 Obelis
 European Authorized Representative Center
 Bis General, Walle, 05, B-1030 Brussels, Belgium