

## INSTRUCTIONS FOR USE

### DEVICE SYSTEM NAME:

#### GREENS SURGICALS (P) LTD. Non-Sterile Bone Plates

##### DEVICE DESCRIPTION:

The bone plates are the devices, used to fasten the bones for the purpose of fixation of fractured bones. The bone plates can be distinguished in terms of their function, i.e. the plate on which this is to be fixed. Generally there are following types of bone plates:

- ▶ Dynamic Compression Type
- ▶ Tubular Type
- ▶ Special (for particular bones)
- ▶ Mini fragment
  
- ▶ The GREENS SURGICALS (P) LTD. Bone Plates are divided mainly into
  - Mini Fragment Plate System
  - Small Fragment Plate System
  - Large Fragment Plate System

The thickness of these plates varies from 1.0mm to 6.0mm  
The width of these plates varies from 3.0mm to 17mm  
The length of these plates varies from 15mm to 360mm  
The number of holes of these plates vary from 2 to 22 holes

The GREENS SURGICALS (P) LTD. Bone Plates are fixed using bone screws as fasteners on the bone, their function, their size and the type of bone they are intended to be used for.

- ▶ Four Types of Screws;
  - Cortical (cortex)
  - Cancellous
  - Malleolar
  - Cannulated

These screw variants can be provided with locking compression thread types as well as regular types.

**Type of Recess :** Hexagonal.  
**Diameter Range :** 1.5mm to 7.3mm.  
**Length Range :** 6mm to 150mm.

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 1 OF 8

## INSTRUCTIONS FOR USE

Adequate knowledge of Surgical Techniques, proper selection and placement of implants and post-operative patient care are important considerations while performing an orthopaedic surgical procedure. Therefore it is essential that the implants are used by a qualified Medical Practitioner only who are surgeons and who know the use of the devices.

**NOTE:** It is the responsibility of the surgeon to discuss, with the patient, the precautions, possible risks, warnings, consequences, complications and adverse reactions which may occur as a result of the surgical procedure and implantation of the device(s). The patient should be informed that the life expectancy of the device is unpredictable and once implanted, successful results cannot be guaranteed if adequate care and precautions are not taken.

**If adverse effects happen, it may necessitate re-operation, revision or removal surgery, arthrodesis or the involved joint and/or amputation of the limb.**

##### MRI COMPATIBILITY:

The GREENS SURGICALS (P) LTD. Bone Plates have not been evaluated for safety and compatibility in the MR environment. The GREENS SURGICALS (P) LTD. Bone Plates have not been tested for heating or migration in the MR environment.

#### I SELECTION OF THE PATIENT

During selection of patients the following factors have to be considered:

- a) **Immunological intolerance**  
Immunological intolerance may occur in some patients. Where material sensitivity is suspected, appropriate foreign body tests should be performed
- b) **Degenerative Diseases**  
In the case of patients suffering from degenerative diseases, this can get aggravated during implantation, and will decrease the expected life of the implant. In such cases surgery can be considered only as a temporary relief.
- c) **Mental illness**  
Mental illness or schizophrenia may cause patients to ignore the limitations and precautions of the implanted material, leading to implants fracture and complication

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 3 OF 8

## INSTRUCTIONS FOR USE

**MATERIALS:** THESE DEVICES CAN BE MADE IN THE FOLLOWING MATERIAL GRADES

**STAINLESS STEEL:** 316L, 316LVM ALLOY

**TITANIUM ALLOY:** TiAl6V4

##### INDICATIONS FOR USE:

GREENS SURGICALS (P) LTD. Bone Plates are provided non-sterile. The bone plates are intended for treating fractures of various bones including the clavicle, pelvis, scapula, long bone (humerus, ulna, radius, femur, tibia and fibula), and small bone (metacarpals, metatarsals, phalanges).

##### CONTRA INDICATIONS:

The implant should not be used in a patient who has currently, or who has a history of:

- ▶ Local or Systemic acute or chronic inflammation.
- ▶ Active infection or inflammation.
- ▶ Suspected or documented metal allergy or intolerance.
- ▶ Symptomatic Arthritis.

##### WARNING:

**Serious Post-operative complications may occur from the use of implant in a patient who-**

- ▶ **Lacks good general physical condition**
- ▶ **Has severe osteoporosis**
- ▶ **Demonstrates anatomical or physiological anomalies.**
- ▶ **Has immunological response, sensitization or hypersensitivity for foreign materials.**
- ▶ **Do Not Re-Use. Do not re-process. The device is not designed for re-processing. Reprocessing of the medical device intended for single use devices may lead to degraded performance or a loss of functionality. The cleaning and disinfection / sterilization for re-processing of these single use devices has not been validated nor any authentic information is available. So re-process of the single use device is not allowed.**

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 2 OF 8

## INSTRUCTIONS FOR USE

#### d) Alcohol and Drug Addiction

The patients, addicted to alcohol and drugs, ignore, during the state of stupor or during the stage of withdrawal, the necessary precautions for the use of implants. This may result in complications of implants fracture.

#### e) Obesity

An obese patient produces abnormal stresses leading to increased load on the implant which ultimately results in failure of the implant.

#### f) Activity

If the patient indulges in an activity involving significant muscular strain in the implanted region, the result may be failure of the implant.

#### II SELECTION OF THE IMPLANT/IMPLANT SYSTEM

For successful implantation, the selection of the proper size, the shape and design of the implant are important. The size and the shape of the human bones puts limitations of and strength of implants during fracture management and re-constructive surgery. The Doctor/Orthopaedic Surgeon has to select the appropriate implant depending upon the application, which should also be free from any apparent corrosion or any manufacturing defect.

#### III HANDLING OF IMPLANT/IMPLANT SYSTEM

The plate should conform to the shape of the bone. Any inequality should be adjusted by accurate shaping of the plate. Hence, plate contouring should be done by the plate benders designed for this purpose.

Care should be taken that there are no scratches or distortions notches, sharp dents or reverse bends at the site of the screw hole. These may cause defects on surface finish and result in improper bonding.

#### IV IMPORTANT

The Plates, even if bent and properly fixed may fatigue and break. Breakage may be due to unstable implant fixation or insufficient support.

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 4 OF 8

## INSTRUCTIONS FOR USE

### V IMPLANT REMOVAL

The surgeon must take the final decision on implant removal. It is recommended that an implant, used as an aid for healing should be removed once its services are over, particularly in younger and more active patients

### VI POST OPERATIVE CARE

A patient must be made aware of the limitations of metallic implants and take precautions to avoid unnecessary stress to implant.

No partial weight bearing or non-weight bearing device can be expected to withstand the unsupported stresses of full weight bearing or excessive muscular activity when there is an un-united fracture. Bone union is a must and the patient must restrict his activities to assist in healing.

### VII SECOND HAND IMPLANTS

Used implants, even if they appear undamaged, may have internal and external Defects. It is possible that individual stress analysis each part fails to reveal the accumulated stress on the metals as a result of use within the body. This may ultimately lead to implant failure.

Every implant must be discarded after use and should never be re-used. It should be bent & then disposed of properly so that it becomes unfit for reuse. While disposing it of, it should be ensured that the discarded implant does not pose any threat to children, stray animals and environment.

### VIII CAUTION

Implant components from one manufacture should not be used with those of another. Implants from each manufacture may have metal and design differences so that the use in conjunction with different devices may lead to inadequate fixation or corrosion of the implant due to generation of peizo currents.

### IX IMPLANT AND PACKAGING DISPOSAL

The packaging material of this device if swallowed, may cause choking Hazards. Therefore, it should be disposed of in such a way that it is out of reach of children and stray animals. Dispose off the implants as per applicable medical practices and as per applicable local, state and country specific regulatory requirements.

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 5 OF 8

## INSTRUCTIONS FOR USE

### INFORMATION

Should any information regarding the implants or their use be required, please contact your representative or distributor or contact the manufacturers direct.

### Details of various Symbols used in Labeling

	Do not reuse
	Batch code
	Date of manufacture YYYY-MM
	Consult Instructions for Use
	Manufacturer
	Use by YYYY-MM
	Non-Sterile
	Catalogue number
	Keep Dry

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 7 OF 8

## INSTRUCTIONS FOR USE

### X STERILITY

**GREENS SURGICALS (P) LTD. Bone Plates are supplied Non-Sterile.**

**Check the integrity of the packaging and labelling before opening the pack.**

**Remove the device from the pack before sterilization.**

**Implants are recommended to be sterilized, using steam autoclaving process regularly used in the hospitals and clinics.**

**The following two methods have been validated and recommended by the company:**

<b>Method:</b>	<b>Steam Sterilization (Autoclaving) Gravity Feed</b>
<b>Temperature</b>	<b>121 Degree Centigrade</b>
<b>Exposure Time:</b>	<b>15 Minutes</b>

### XI USE

Qualified & specialized trained surgeons only should use this device.

### STORAGE

Store in a dry place.

### PRODUCT INFORMATION DISCLOSURE:

GREENS SURGICALS (P) LTD. HAS EXERCISED REASONABLE CARE IN THE SELECTION OF MATERIALS AND MANUFACTURE OF THESE PRODUCTS. GREENS SURGICALS (P) LTD. EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. GREENS SURGICALS (P) LTD. WILL NOT BE LIABLE FOR ANY INCIDENT OR CONSEQUENTIAL LOSS, DAMAGE OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM THE USE OF THE DEVICES. GREENS SURGICALS (P) LTD. NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE DEVICES.

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 6 OF 8

## INSTRUCTIONS FOR USE

	Keep Away from Sunlight
	Do not use if package is damaged
	EC Representative Address
	Material Grade Stainless Steel like 316L, 316LVM, Titanium TiA6V4
	Quantity in the package

Manufactured By:  
**GREENS SURGICALS (P) LTD.**  
FDA, ISO 13485/2003 Certified  
Reg No: 516-612, C/D, Suli Industrial Estate,  
Mumbai, Maharashtra 401115 (Goregaon)  
Phone: +91-02267-294569, 294569, Toll Free: +91-02267-294880  
Web Site: www.greensurgicals.com  
Email: info@greensurgicals.com

**Obelis**  
European Authorized Representative Center  
Bd General Waha, 25, B-1030 Brussels, Belgium

SECTION : GREENS SURGICALS (P) LTD. /IFU/BONE PLATES  
REV.01 Dated: 17-07-2015 PAGE 8 OF 8