



(VETERINARY)

PLATES

(2014-15)



GREENS SURGICALS (P) LTD

Plot No. 508-512, G.I.D.C. Savli Industrial Estate, Manjusar, Vadodara - 391775, (Gujarat), INDIA

TEL : +91-2667 - 264888 & +91-265-2783892

EMAIL: info@greensurgicals.com

WEB : www.greensurgicals.com


1 LC-DCP PLATES

2.0mm/1.5mm LC-DCP Plates

Width: 5.5mm

Distance between holes 6.0mm

For use with 1.5mm or 2.0mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.2 mm
500.201	4	25	
500.202	5	31	
500.203	6	37	
500.204	7	43	

Part #	Holes	Length mm	Thickness - 1.5 mm
500.211	6	37	
500.212	7	43	
500.213	8	49	
500.214	9	55	
500.215	10	61	
500.216	12	73	
500.217	14	85	

2.4mm/2.0mm LC-DCP Plates

Width: 6.5mm

Distance between holes 8.0mm LC-DCP PLATE

For use with 2.0mm or 2.4mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 2.0 mm
500.221	8	65	
500.222	9	73	
500.223	10	81	
500.224	12	97	
500.225	14	113	


Part #	Holes	Length mm	Thickness - 1.7 mm
500.231	4	33	
500.232	5	41	
500.233	6	49	
500.234	7	57	
500.235	8	65	

2.7mm LC-DCP Plates

Width: 7.5mm

Distance between holes 9.0mm

For use with 2.7mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 2.6 mm
500.241	4	37	
500.242	5	46	
500.243	6	55	
500.244	7	64	
500.245	8	73	
500.246	9	82	
500.247	10	91	
500.248	11	100	
500.249	12	109	
500.250	14	127	
500.251	16	145	

3.5mm LC-DCP Plates

Width: 11.0mm

Distance between holes 13.0mm

For use with 3.5mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 3.3 mm
500.261	2	28	
500.262	3	41	
500.263	4	54	
500.264	5	67	
500.265	6	80	
500.266	7	93	
500.267	8	106	
500.268	9	119	
500.269	10	132	
500.270	11	145	
500.271	12	158	
500.272	14	184	
500.273	15	197	
500.274	16	210	
500.275	18	236	
500.276	20	262	
500.277	22	288	

3.5mm Broad LC-DCP Plates

Width: 13.5mm

Distance between holes 13.0mm

For use with 3.5mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 4.2 mm
500.281	7	94	
500.282	8	107	
500.283	9	120	
500.284	10	133	
500.285	11	146	
500.286	12	159	
500.287	13	172	
500.288	14	185	
500.289	15	198	
500.290	16	211	
500.291	17	224	
500.292	18	237	
500.293	20	263	
500.294	22	289	

4.5mm Narrow LC-DCP Plates

Width: 13.5mm

Distance between holes 18.0mm

For use with 4.5mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 5.2 mm
500.301	2	34	
500.302	3	52	
500.303	4	70	
500.304	5	88	
500.305	6	106	
500.306	7	124	
500.307	8	142	
500.308	9	160	
500.309	10	178	
500.310	11	196	
500.311	12	214	
500.312	13	232	
500.313	14	250	
500.314	15	268	
500.315	16	286	
500.316	18	323	
500.317	20	358	
500.318	22	394	

4.5mm Broad LC-DCP Plates

Width: 17.5mm

Distance between holes 18.0mm

For use with 4.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 5.2 mm
500.325	6	106	
500.326	7	124	
500.327	8	142	
500.328	9	160	
500.329	10	178	
500.330	11	196	
500.331	12	214	
500.332	14	250	
500.333	16	286	
500.334	18	322	
500.335	20	358	
500.336	22	394	

2 LCP PLATES


2.0mm/1.5mm LCP Plates


Width: 5.5mm

Distance between centre of holes 7.0mm

Distance between central holes 4.6mm

For use with 2.0mm Locking, 1.5mm Cortex, and 2.0mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.2 mm
500.341	4	27	
500.342	5	34	
500.343	6	41	
500.344	7	48	

Part #	Holes	Length mm	Thickness - 1.5 mm
500.351	6	41	
500.352	7	48	
500.353	8	55	
500.354	9	62	
500.355	10	69	
500.356	12	83	
500.357	14	97	


2.4mm LCP Plates


Width: 6.5mm

Distance between centre of holes 9.0mm

Distance between central holes 5.6mm

For use with 2.4mm Locking, 2.0mm Cortex, and 2.4mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.7 mm
500.361	4	32	
500.362	5	40	
500.363	6	48	
500.364	7	56	
500.365	8	64	

Part #	Holes	Length mm	Thickness - 2.0 mm
500.371	8	64	
500.372	9	72	
500.373	10	80	
500.374	12	96	
500.375	14	112	


2.7mm LCP Plates

Width: 7.5mm

Distance between centre of holes 9.0mm

Distance between central holes 6.0mm

For use with 2.7mm Locking, and 2.7mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 2.6 mm
500.381	4	36	
500.382	5	45	
500.383	6	54	
500.384	7	63	
500.385	8	72	
500.386	9	81	
500.387	10	90	
500.388	12	108	
500.389	14	126	
500.390	16	144	


3.5mm LCP Plates

Width: 11.0mm

Distance between centre of holes 13.0mm

Distance between central holes 9.4mm

For use with 3.5mm Locking, 3.5mm Cortex, and 4.0mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 3.3 mm
500.395	2	27	
500.396	3	40	
500.397	4	53	
500.398	5	66	
500.399	6	79	
500.400	7	92	
500.401	8	105	
500.402	9	118	
500.403	10	131	
500.404	11	144	
500.405	12	157	
500.406	13	170	
500.407	14	183	
500.408	15	196	
500.409	16	209	
500.410	18	235	
500.411	20	261	
500.412	22	287	


3.5mm Broad LCP Plates

Width: 13.5mm

Distance between centre of holes 13.0mm

Distance between central holes 9.4mm

For use with 3.5mm Locking, 3.5mm Cortex, and 4.0mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 4.2 mm
500.421	7	94	
500.422	8	107	
500.423	9	120	
500.424	10	133	
500.425	11	146	
500.426	12	159	
500.427	13	172	
500.428	14	185	
500.429	15	198	
500.430	16	211	
500.431	17	224	
500.432	18	237	
500.433	20	263	
500.434	22	289	


4.5mm Narrow LCP Plates

Width: 13.5mm

Distance between centre of holes 18.0mm

Distance between central holes 12.8mm

For use with 4.0mm and 5.0mm Locking, 4.5mm Cortex, and 6.5mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 4.2 mm
500.441	3*	66*	
500.442	4	71	
500.443	5	89	
500.444	6	107	
500.445	7	125	
500.446	8	143	
500.447	9	161	
500.448	10	179	
500.449	11	197	
500.450	12	215	
500.451	14	251	
500.452	16	287	

* for proximal interphalangeal arthrodesis.


4.5mm Broad LCP Plates

Width: 17.5mm

Distance between centre of holes 18.0mm

Distance between central holes 12.8mm

For use with 4.0mm and 5.0mm Locking, 4.5mm Cortex, and 6.5mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 5.2 mm
500.461	6	107	
500.462	7	125	
500.463	8	143	
500.464	9	161	
500.465	10	179	
500.466	11	197	
500.467	12	215	
500.468	13	233	
500.469	14	251	
500.470	15	269	
500.471	16	287	
500.472	18	323	


5.5mm Broad LCP Plates

Width: 17.5mm

Distance between centre of holes 18.0mm

Distance between central holes 12.8mm

For use with 4.0mm and 5.0mm Locking, 4.5mm and 5.5mm Cortex, and 6.5mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 6.0 mm
500.481	10	180	
500.482	11	198	
500.483	12	216	
500.484	13	234	
500.485	14	252	
500.486	15	270	
500.487	16	288	
500.488	17	306	
500.489	18	324	


3.5mm LCP Recon Plates

Width: 10.1mm

Distance between centre of holes 14.0mm

Distance between central holes 10.0mm

For use with 3.5mm Locking, 3.5mm Cortex, and 4.0mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 2.8 mm
500.495	4	56	
500.496	5	70	
500.497	6	84	
500.498	7	98	
500.499	8	112	
500.500	9	126	
500.501	10	140	
500.502	11	154	
500.503	12	168	

3.5mm LCP T-Plates 3 Holes Head

Width: 10.0mm

Distance between holes, 1.5mm head holes accept 3.5mm locking, 2.7mm and 3.5mm Cortex Screws
Shaft holes accept 3.5mm Locking, 3.5mm Cortex, and 4.0mm Cancellous Bone Screws.


Right angle

Part #	Holes	Length mm	Thickness - 2.8 mm
500.511	3	50	
500.512	4	57	
500.513	5	67	
500.514	6	77	
500.515	7	87	
500.516	8	97	

Notched Head T-Plates 2 Holes Head


2.0mm LCP Notched Head T-Plate, 7 holes shaft, 54mm

For use with 2.0mm Locking and 2.0mm Cortex Screws.

Part #	Description	Thickness - 1.3 mm
500.521	Width of Head: 10.0mm Width of Shaft: 5.0mm	

2.4mm LCP Notched Head T-Plate, 7 holes shaft, 59mm

For use with 2.4mm Locking and 2.4mm Cortex Screws.

Part #	Description	Thickness - 1.3 mm
500.526	Width of Head: 12.0mm Width of Shaft: 6.4mm	

2.7mm LCP Notched Head T-Plate, 7 holes shaft, 66mm

For use with 2.7mm Locking and 2.7mm Cortex Screws.

Part #	Description	Thickness - 1.5 mm
500.531	Width of Head: 14.2mm Width of Shaft: 7.0mm	

2.0mm Ti Mini Locking Recon Plate

Width: 4.8mm

Hole spacing, 6.5mm

For use with 2.0mm Ti Locking and 2.0mm Ti Cortex Screws.

Part #	Holes	Type	Thickness - 1.0 mm
500.544	20	Adaption	Picture coming soon
500.536	4	Narrow	
500.537	6	Narrow	

2.0mm Ti Intermediate Locking Recon Plate

Width: 5.0mm

Hole spacing, 6.5mm

For use with 2.0mm Ti Locking and 2.0mm Ti Cortex Screws.

Part #	Holes	Thickness - 1.3 mm
500.553	6	Picture coming soon
500.556	12	

2.0mm Ti Large Locking Recon Plate

Width: 6.5mm

Hole spacing, 8.0mm

For use with 2.0mm Ti Locking and 2.0mm Ti Cortex Screws.

Part #	Holes	Thickness - 1.5 mm
500.565	3x3	Picture coming soon
500.567	6	
500.570	12	
500.574	20	

2.0mm Ti Extra-Large Locking Recon Plate

Width: 6.5mm

For use with 2.0mm Ti Locking and 2.0mm Ti Cortex Screws.

Part #	Holes	Thickness - 2.0 mm
500.585	12	Picture coming soon
500.588	20	



2.4mm Ti Locking Recon Plate

Width: 8.0mm

Hole spacing, 8.0mm

For use with 2.4mm and 3.0mm Ti Locking and 2.4mm Ti Cortex Screws.

Part #	Holes	Thickness - 2.5 mm
500.591	14	Picture coming soon
500.595	24	

3 DCP PLATES


1.5mm/2.0mm DCP Plates


Width: 5.0mm

Distance between centre of holes 6.0mm

Distance between central holes 7.0mm

For use with 1.5mm and 2.0mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.5 mm
500.601	6	32	
500.602	7	37	
500.603	8	42	

Part #	Holes	Length mm	Thickness - 1.0 mm
500.611	4	22	
500.612	5	27	
500.613	6	32	


2.7mm DCP Plates


Width: 8.0mm

Distance between centre of holes 8.0mm

Distance between central holes 12.0mm

For use with 2.7mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 2.0 mm
500.621	2	20	
500.622	4	36	
500.623	5	44	
500.624	6	52	

Part #	Holes	Length mm	Thickness - 2.5 mm
500.631	7	60	
500.632	8	68	
500.633	9	76	
500.634	10	84	
500.635	12	100	


3.5mm DCP Plates

Width: 10.0mm

Distance between centre of holes 12.0mm

Distance between central holes 16.0mm

For use with 3.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 3.3 mm
500.641	2	25	
500.642	3	37	
500.643	4	49	
500.644	5	61	
500.645	6	73	
500.646	7	85	
500.647	8	97	
500.648	9	109	
500.649	10	121	
500.650	12	145	
500.651	14	169	
500.652	16	194	
500.653	18	217	
500.654	20	242	
500.655	22	265	


3.5mm Broad DCP Plates

Width: 12.0mm

Distance between centre of holes 12.0mm

Distance between central holes 16.0mm

For use with 3.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 3.6 mm
500.661	7	86	
500.662	8	98	
500.663	9	110	
500.664	10	122	
500.665	11	134	
500.666	12	146	
500.667	13	158	
500.668	14	170	
500.669	16	194	
500.670	18	218	
500.671	20	242	
500.672	22	266	


4.5mm Narrow DCP Plates

Width: 12.0mm

Distance between centre of holes 16.0mm

Distance between central holes 25.0mm

For use with 4.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 3.8 mm
500.681	2	39	
500.682	3	55	
500.683	4	71	
500.684	5	87	
500.685	6	103	
500.686	7	119	
500.687	8	135	
500.688	9	151	
500.689	10	167	
500.690	11	183	
500.691	12	199	
500.692	13	215	
500.693	14	231	
500.694	15	247	
500.695	16	263	
500.696	18	295	
500.697	20	326	
500.698	22	358	
500.699	24	390	


4.5mm Broad DCP Plates

Width: 16.0mm

Distance between centre of holes 16.0mm

Distance between central holes 25.0mm

For use with 4.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 4.8 mm
500.705	6	103	
500.706	7	119	
500.707	8	135	
500.708	9	151	
500.709	10	167	
500.710	12	199	
500.711	14	231	
500.712	16	263	
500.713	18	295	
500.714	20	326	
500.715	22	358	
500.716	24	390	
500.717	26	423	

2.0mm Acetabular Plate

Width: 5.0mm

Hole spacing 15° (5mm)

For use with 2.0mm Cortex Screws.

Part #	Holes	Length mm	
500.721	6	25	

2.7mm Acetabular Plate

Width: 8.0mm

Hole spacing 20.5° (9.6mm)

For use with 2.7mm Cortex Screws.

Part #	Holes	Length mm	
500.725	6	45	


2.7mm Triple Osteotomy Plates, 45°

Width: 8.0mm

Distance between holes 8mm

Distance between central holes 13.0mm

For use with 2.7mm Cortex Screws.

Part #	Holes	Length mm	
500.731	5	45 (left)	
500.732	6	53 (left)	
500.733	5	45 (right)	
500.734	6	53 (right)	


3.5mm Triple Osteotomy Plates, 45°

Width: 10.0mm

Distance between holes 16.0mm

Distance between central holes 17.0mm

For use with 3.5mm Cortex Screws.


Part #	Holes	Length mm	
500.741	5	63 (left)	
500.742	6	75 (left)	
500.743	5	63 (right)	
500.744	6	75 (right)	

4.0mm/3.5mm Double Hook Plate

Width: 12.0mm

Distance between holes 16.0mm

For use with 3.5mm Cortex and 4.0mm Cancellous Bone Screws.

Part #	Holes	Length mm	Thickness - 3.8 mm
500.751	5	84	
500.753	7	116	

4 ROUND HOLE PLATES

1.5mm Straight Plate

Width: 3.8mm

Distance between centre of holes 5.0mm

For use with 1.5mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 0.8 mm
500.762	20	100	

2.0mm Straight Plate

Width: 5.0mm

Distance between centre of holes 6.0mm

For use with 2.0mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.0 mm
500.765	3	17	
500.766	4	23	
500.767	5	29	
500.768	6	35	
500.773	20	100	

2.0mm/1.5mm Cut-To-Length Plate

Width: 7.0mm

Distance between centre of holes 6.0mm

For use with 1.5mm and 2.0mm Cortex Screws

May be cut to length to accommodate a variety of fracture situations with Plate Cutting Pliers.

Part #	Holes	Length mm	Thickness - 1.0 mm
500.776	50	300	

2.7mm/2.0mm Cut-To-Length Plate

Width: 7.0mm

Distance between centre of holes 6.0mm

For use with 2.0mm and 2.7mm Cortex Screws

May be cut to length to accommodate a variety of fracture situations with Plate Cutting Pliers.


Part #	Holes	Length mm	Thickness - 1.5 mm
500.778	50	300	

1.5mm T-Plates, 9 Holes in Shaft, 50mm

Width: 3.8mm

Distance between centre of holes 5.0mm

For use with 1.5mm Cortex Screws.

Part #	Head Holes	Holes	Thickness - 0.9 mm
500.781	3	9	
500.782	4	9	

2.0mm T-Plate 2 Holes, 18mm

Width: 5.0mm

Distance between centre of holes 6.0mm

For use with 2.0mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 1.0 mm
500.786	2	18	

2.0mm T-Plates, 9 Holes in Shaft, 50mm

Width: 5.0mm


Distance between centre of holes 5.0mm

For use with 2.0mm Cortex Screws.

Part #	Head Holes	Holes	Thickness - 1.0 mm
500.790	3	9	
500.791	4	9	


2.0mm/1.5mm T-Plates, 2 Holes in Head

Width: 5.0mm
 Width of T-Head, 9.0mm
 Width of plate shaft, 5.0mm
 For use with 2.0mm and 1.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.5 mm
500.795	4	34	
500.799	8	62	

2.7mm T-Plate, 2 Holes in Head

Width: 7.0mm
 Distance between centre of holes 8.0mm
 For use with 2.7mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.2 mm
500.800	3	32	

3.5mm/2.7mm T-Plates, 2 Holes in Head


Width: 5.0mm
 Width of T-Head, 18.0mm
 Width of plate shaft, 9.0mm
 For use with 2.7mm and 3.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 2.0 mm
500.803	5	65	
500.805	7	85	
500.807	9	105	

3.5mm T-Plates, right angle

Width: 10.0mm

For use with 3.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.2 mm
500.810	3*	50	
500.811	5*	67	
500.812	4**	56	
500.813	6**	78	
500.814	7**	90	
500.815	8**	101	
500.816	9**	112	


* 3holes head

** 4holes head

T-Plates, 2 holes in head

Width of shaft: 17.3mm

For use with 4.5mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 2.0 mm
500.821	3	68	
500.822	4	84	
500.823	5	100	
500.824	6	116	
500.825	8	148	
500.826	11	198	
500.827	12	210	

2.0mm L-Plates, 2 Shaft holes

Width: 5.0mm

Distance between holes, 6.0mm

For use with 2.0mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 1.0 mm
500.831	2	17 left	
500.832	2	17 right	
500.833	2	20 left oblique	
500.834	2	20 right oblique	

2.7mm L-Plates, 3 Shaft holes

Width: 7.0mm

Distance between centre of holes, 8.0mm

For use with 2.7mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.0 mm
500.838	3	32 left	
500.839	3	32 right	
500.840	3	34 left oblique	
500.841	3	34 right oblique	

1.5mm H-Plate

Length: 9.0mm

Width: 8.2mm

Distance between centre of holes, 5.0mm length and 4.2mm width

For use with 1.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 1.0 mm
500.845	4	4.2	

2.0mm H-Plate

Length: 12.0mm

Width: 11.0mm

Distance between centre of holes, 7.0mm length and 6.0mm width

For use with 2.0mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 1.0 mm
500.848	4	6.0	

2.7mm Reconstruction Plates

Width: 8.0mm

Distance between centre of holes, 8.0mm

For use with 2.7mm Cortex Screws.


Part #	Holes	Length mm	Thickness - 3.3 mm
500.851	5	40	
500.852	6	48	
500.853	8	64	
500.854	10	80	
500.855	12	96	
500.856	14	112	
500.857	16	128	
500.858	18	144	
500.859	20	160	
500.860	22	176	
500.861	24	192	

3.5mm Reconstruction Plates

Width: 10.0mm

Distance between centre of holes, 12.0mm

For use with 3.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 2.8 mm
500.865	4	46	
500.866	5	58	
500.867	6	70	
500.868	7	82	
500.869	8	94	
500.870	9	106	
500.871	10	118	
500.872	11	130	
500.873	12	142	
500.874	13	154	
500.875	14	166	
500.876	15	178	
500.877	16	190	
500.878	18	214	
500.879	20	238	
500.880	22	263	


3.5mm Curved Reconstruction Plates

Width: 10.0mm

Distance between centre of holes, 12.0mm

Radius of curvature, 100mm

For use with 3.5mm Cortex Screws.

Part #	Holes	Length mm	Thickness - 3.6 mm
500.885	6	70	
500.886	8	94	
500.887	10	118	
500.888	12	142	
500.889	14	166	
500.890	16	190	
500.891	18	214	


4.5mm Narrow Lengthening Plates

Width: 12.0mm

Distance between centre of holes, 12.0mm

Radius of curvature, 100mm

For use with 4.5mm Cortex Screws, 6.5mm. Cancellous Bone Screws can be used in the end holes.

Part #	Holes	Lengthens	Thickness - 3.8 mm
500.895	8	30 mm 135 mm	
500.896	8	40 mm 145 mm	
500.897	8	50 mm 155 mm	
500.898	8	60 mm 165 mm	
500.899	8	70 mm 175 mm	
500.900	8	80 mm 185 mm	


4.5mm Broad Lengthening Plates

Width: 16.0mm

Distance between centre of holes, 12.0mm

Radius of curvature, 100mm

For use with 4.5mm Cortex Screws, 6.5mm. Cancellous Bone Screws can be used in the end holes.

Part #	Holes	Lengthens	Thickness - 4.8 mm
500.905	8	30 mm 135 mm	
500.906	8	40 mm 145 mm	
500.907	8	50 mm 155 mm	
500.908	8	60 mm 165 mm	
500.909	10	50 mm 179 mm	
500.910	10	60 mm 189 mm	
500.911	10	70 mm 199 mm	
500.912	10	80 mm 209 mm	
500.913	10	90 mm 219 mm	
500.914	10	100 mm 229 mm	
500.915	10	110 mm 239 mm	
500.916	10	120 mm 249 mm	

5 TPLO PLATES

2.0mm TPLO Plates

Width: 5.5mm

Head holes accept 2.0mm Locking and 2.0mm Cortex Screws

Shaft holes accept 2.0mm Cortex Screws

Shaft combination hole accepts 2.0mm Locking or 2.0mm Cortex Screws.

Part #	Holes	Side	Thickness - 2.3 mm
500.921	3	left	
500.922	3	right	

2.4mm TPLO Plates

Width: 6.5mm

Head holes accept 2.4mm Locking and 2.4mm Cortex Screws

Shaft holes accept 2.4mm Cortex Screws

Shaft combination hole accepts 2.4mm Locking or 2.4mm Cortex Screws.

Part #	Holes	Side	Thickness - 2.3 mm
500.926	3	left	
500.927	3	right	


2.7mm TPLO Plates

Width: 8.0mm

Head holes accept 2.7mm Locking and 2.7mm Cortex Screws

Shaft holes accept 2.7mm Cortex Screws

Shaft combination hole accepts 2.7mm Locking or 2.7mm Cortex Screws.

Part #	Holes	Side	Thickness - 2.5 mm
500.931	3	left	
500.932	3	right	


3.5mm TPLO Plates

Width: 11.0mm

Head holes accept 3.5mm Locking and 4.0mm Cortex Screws

Shaft holes accept 3.5mm Cortex Screws

Shaft combination hole accepts 3.5mm Locking or 3.5mm Cortex Screws.

Part #	Holes	Side	Thickness - 3.7 mm
500.936	3	left	
500.937	3	right	

3.5mm Broad TPLO Plates



Width: 13.5mm

Head holes accept 3.5mm Locking, 3.5mm Cortex and 4.0mm Cancellous Bone Screws

Shaft holes accept 3.5mm Cortex Screws

Combination holes accept 3.5mm Locking or 3.5mm Cortex Screws.

Part #	Holes	Side	Thickness - 3.7 mm
500.941	4	left	
500.942	4	right	

3.5mm Small TPLO Plates

Width: 11.0mm

Head holes accept 3.5mm Locking, 3.5mm Cortex and 4.0mm Cancellous Bone Screws

Shaft holes accept 3.5mm Cortex Screws

Combination holes accept 3.5mm Locking or 3.5mm Cortex Screws.

Part #	Holes	Side	Thickness - 3.7 mm
500.946	3	left	
500.947	3	right	